

Dialysis

GeBAflex-tube Kit

Mini GeBAflex-tube Vol. 10-250 μ l

Product	Units	MWCO	Catalog
Mini GeBAflex-tube	10	6000 - 8000	D070-6-10
Mini GeBAflex-tube	30	6000- 8000	D070-6-30
Mini GeBAflex-tube	50	6000-8000	D070-6-50
Mini GeBAflex-tube	100	6000-8000	D070-6-100
Mini GeBAflex-tube	10	12000-14000	D070-12-10
Mini GeBAflex-tube	30	12000-14000	D070-12-30
Mini GeBAflex-tube	50	12000-14000	D070-12-50
Mini GeBAflex-tube	100	12000-14000	D070-12-100
Mini GeBAflex-tube	5	25,000	D070-25-5

Midi GeBAflex-tube Vol. 50-800 μ l

Product	Units	MWCO	Catalog
Midi GeBAflex-tube	5	1000	D015-1-5
Midi GeBAflex-tube	10	3500	DO10
Midi GeBAflex-tube	30	3500	DO12
Midi GeBAflex-tube	50	3500	D010-50
Midi GeBAflex-tube	100	3500	D010-100
Midi GeBAflex-tube	10	6000-8000	DO20
Midi GeBAflex-tube	30	6000-8000	DO22
Midi GeBAflex-tube	50	6000-8000	D020-50
Midi GeBAflex-tube	100	6000-8000	D020-100

Maxi GeBAflex-tube Vol. 0.1-3ml

Product	Units	MWCO	Catalog
Maxi GeBAflex-tube	5	3500	D030
Maxi GeBAflex-tube	15	3500	D035
Maxi GeBAflex-tube	50	3500	D030-50
Maxi GeBAflex-tube	100	3500	D030-100
Maxi GeBAflex-tube	5	6000-8000	DO40
Maxi GeBAflex-tube	15	6000-8000	DO45
Maxi GeBAflex-tube	50	6000-8000	D040-50
Maxi GeBAflex-tube	100	6000-8000	D040-100
Maxi GeBAflex-tube	5	12000-14000	DO50
Maxi GeBAflex-tube	15	12000-14000	DO55
Maxi GeBAflex-tube	50	12000-14000	D050-50
Maxi GeBAflex-tube	100	12000-14000	D050-100
Maxi GeBAflex-tube	5	25,000 MWCO	D060-25-5
Maxi GeBAflex-tube	5	50,000 MWCO	D060-50-5



Mini GeBAflex-tube & float



Midi GeBAflex-tube & float



Maxi GeBAflex-tube & float.
with 2 and 3 ml cap

<p>Shelf-life</p> <p>12 month</p>	<p>On line Support</p> <p>geba@geba.org</p>	<p>Storage</p> <p>15-25^oC</p>	<p>Storage</p> <p>1000 MwCO 25000 MwCO 50000 MwCO</p> <p>+4^oC</p>
-----------------------------------	---	--	--

Simply Better

Distributed by:


CliniSciences Group

Dialysis

GeBAflex-tube Kit

Mega GeBAflex-tube Vol. range 3-20ml

Product	Units	MWCO	Catalog
Mega GeBAflex-tube 10ml cap	5	1000	D080-1/10-5
Mega GeBAflex-tube 15 ml cap	5	1000	D080-1/15-5
Mega GeBAflex-tube 20 ml cap	5	1000	D080-1/20-5
Mega GeBAflex-tube 20 ml cap	10	3500	DO80-3,5/20-10
Mega GeBAflex-tube 15 ml cap	10	3500	DO80-3,5/15-10
Mega GeBAflex-tube 10 ml cap	10	3500	DO80-3,5/10-10
Mega GeBAflex-tube 20 ml cap	10	6000-8000	DO80-6/20-10
Mega GeBAflex-tube 15 ml cap	10	6000-8000	DO80-6/15-10
Mega GeBAflex-tube 10 ml cap	10	6000-8000	DO80-6/10-10
Mega GeBAflex-tube 20 ml cap	10	12000-14000	DO80-12/20-10
Mega GeBAflex-tube 15 vml cap	10	12000-14000	DO80-12/15-10
Mega GeBAflex-tube 10 ml cap	10	12000-14000	DO80-12/10-10



Mega GeBAflex-tube & float with 10, 15, 20 ml cap

- ◆ **GeBAflex-tube** takes dialysis & buffer exchange process a step further in human engineering, in high throughput efficiency for volume ranging from 10 μ l to 20 ml.
- ◆ **GeBAflex-tube** is ready to use, membrane is RNAs, DNAs, sulfur & heavy metal free and EDTA treated.
- ◆ Wide membrane range from 1000 MWCO to 50 000 MWCO available
- ◆ Standard pipette tip is used for loading & unloading GeBAflex-tube.
- ◆ Volume Recovery is typically in the 98% range.
- ◆ Single float T006-63 can process upto 63 Mini GeBAflex-tube in a single beaker
Single float T002-35 can process upto 35 Midi GeBAflex-tube in a single beaker
Single float T003-20 can process upto 20 Maxi GeBAflex-tube in a single beaker
- ◆ GeBAflex-tube caps option improve dialysis effect by maximizing liquid surface membrane contact
- ◆ GeBAflex-tube kit contain handbook instruction & float

High through put dialysis & buffer exchange



Single float T002-35 can process upto 35 Midi GeBAflex-tube in a single beaker



Simply Better

Distributed by:

CliniSciences Group

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polaska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences Limited
Address: Marktgasse 18 8302 Kloten -
Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>